

## Instructions for collecting single cells (nuclei)

Version: TISP161101

Single cell collection can be performed using multiple methods, for example, mouth pipette [1], Fluorescence-Activated Cell Sorting (FACS) [2] or CellRaft system [3].

### Sample requirement

Single cells (nuclei) in individual 0.1 or 0.2-ml PCR tubes with 2.5  $\mu$ l 1 x PBS (excluding calcium and magnesium).

Store at -80°C.

### Instructions for single cell (nucleus) isolation using FACS

#### *Before FACS*

1. Prepare PCR tube strips (with cap strips) with 2.5  $\mu$ l PBS.  
Caution: Avoid any contamination.  
Suggestion: In PCR workstation or biological safety cabinets.
2. Quick-spin tubes to bring the PBS to their bottom.

#### *After FACS*

1. Deposit single cells (nuclei) into PCR strips.  
Caution: Ensure that cells are placed in the middle of the tubes.
2. Quick-spin tubes.
3. Freeze immediately on dry ice.
4. Store at -80°C.

For further assistance, please contact [support@singulomics.com](mailto:support@singulomics.com).

### References

1. Spits, C., et al., *Whole-genome multiple displacement amplification from single cells*. Nat Protoc, 2006. **1**(4): p. 1965-70.
2. Basu, S., et al., *Purification of specific cell population by fluorescence activated cell sorting (FACS)*. J Vis Exp, 2010(41).
3. <http://cellmicrosystems.com/>.